

MS-Compatible Silver staining of proteins in SDS- and native-PAGE

A newer Gel staining protocol

Chevallet M et al. (2008) Sweet silver: A formaldehyde-free silver staining using aldoses as developing agents, with enhanced compatibility with mass spectrometry. *Proteomics*, 8:4853

An older Gel staining protocol (adapted from Blum et al. *Electrophoresis*, 8, 93-99, 1987)

Prepare all solutions fresh!

This is an easy staining protocol and is compatible with mass spectrometric protein analysis. The stain is useful for protein concentrations ranging from <1ug to >1ng. If the staining does not work, for whatever reason, then gels can easily be destained and restained again (see below).

| | solution | | V = 100 ml | operation | time |
|----|---------------|--|---------------------------------------|---|------------|
| 1 | sol A | 50% methanol, 10% acetic acid | | fix | 30 min |
| 2 | sol B | 50% methanol | | incubate | 15 min |
| 3 | | milli-Q H ₂ O | | Wash, several changes of H ₂ O | 1 hour |
| 4 | sol. C | Sodium thiosulphate (Na₂S₂O₃·5H₂O) | 0.2 g/L fresh! | Incubate | 120 sec |
| 5 | | milli-Q H ₂ O | | wash 3 times | 2 x 30 sec |
| 6 | sol D | silver nitrate (AgNO ₃) | 0.2 g/100 mL | Incubate | >25 min |
| 7 | | milli-Q H ₂ O | | wash 3 times | 3 x 60 sec |
| 8 | sol E | Sodium carbonate (Na ₂ CO ₃) 37% HCOH Na ₂ S ₂ O ₃ x 5H ₂ O (Sol C) | 3 g/100 mL 25ul/100mL 2mL/100mL | develop | max 10 min |
| 9 | sol F | Na ₂ -EDTA | (14g / L) | stop develop | 10 min |
| 10 | | milli-Q H ₂ O | | wash | |

Destaining protocol

- Dissolve 0.4g Potassium ferricyanide (K₃Fe(CN)₆) in 200ml Sodium thiosulphate (0.2 g/L, **sol. C** above)

Destain:

- destain in ferricyanide solution (**sol. C**) with several changes until no bands are visible; the gel will have a yellow hue!
- wash gel 4-5 times for 15 min with milli-Q H₂O until gel is transparent and has no background colour.
- stain gel again starting with **sol. C** from gel staining protocol.